

# First Measurement of Brain-Derived Tau in Plasma Using a Novel Luminescence Immunoassay.

Steffen Busch<sup>1</sup>, Ingolf Lachmann<sup>2</sup>, Oliver Schmidt<sup>3</sup>, Hermann Esselmann<sup>4</sup>, Barbara Morgado<sup>4</sup>, Oliver Schmidt<sup>3</sup>, Jens Wiltfang<sup>4</sup>, Peter Findeisen<sup>1</sup>

<sup>1</sup>MVZ Labor Dr. Limbach & Kollegen eGbR, Im Breitspiel 16, 69126 Heidelberg, Germany; <sup>2</sup>Roboscreen GmbH, Hohmannstraße 7, 04129 Leipzig, Germany; <sup>3</sup>Tecan, IBL International GmbH, Flughafenstraße 52A, 22335 Hamburg, Germany, <sup>4</sup>Psychiatry of University Medical Center Göttingen, Von-Siebold-Str. 5, 37075 Göttingen, Germany,

## Background

Early detection of Alzheimer's disease (AD) risk requires sensitive and specific biomarkers measurable in blood. We evaluated a new Luminescence Immunoassay (LUM) for the quantification of brain-derived tau (BD-Tau) in plasma or serum. The method includes an upstream enrichment step, termed Neuro-IP, employing a monoclonal antibody specific for BD-Tau.




## Methods

Three panels of plasma samples obtained from the University of Göttingen, MVZ Labor Limbach, and a commercial source, all classified as characterized by CSF levels of A $\beta$ 40/42, p181Tau, and hTau, measured with the Lumipulse system, were used. The subjects were classified according to the clinical cutoff point into the diagnostic groups A $\beta$ -positive (AD, n=46) and A $\beta$ -negative (CTRL, n=22). BD-Tau was extracted from a 200  $\mu$ l sample using BD-specific Neuro-IP according to the protocol. The eluates were quantified using BD-Tau-specific Luminescence Immunoassay.

## Conclusion

For the first time, our laboratory investigated the measurement of BD-Tau on routine plasma samples using a standard luminescence microplate immunoassay after Neuro-IP enrichment. This immunoprecipitation step minimizes matrix effects and enables sensitive BD-Tau detection. The BD-Tau LUM assay provides information that complements existing plasma tau biomarkers and may improve the interpretation of patient status in the context of Alzheimer's disease risk assessment.

## Results



(A) BD-specific Neuro-IP with subsequent measurement of BD-Tau concentration using the BD-Tau LUM assay showed significant differentiation between the AD group and control group ( $p < 0.001$ ). (B) The ROC analysis for differentiation of the two diagnostic groups showed an AUC of 0.922 for the BD-Tau LUM assay ( $p < 0.001$ ). (C) Duplicate testing of BD-Tau in plasma from n=51 patients shows perfect positive correlation for the two runs ( $r = 1.00$ ;  $p < 0.001$ ). Furthermore, plasma pTau217 levels showed a strong correlation with BD-Tau concentrations determined by the LUM assay (not shown). Notably, BD-Tau demonstrated the strongest association with CSF A $\beta$  pathology among all measured parameters.