

SUMMERTIME IS TICK TIME!

TBE/FSME IGG PLUS ELISA & TBE/FSME IGM ELISA FROM TECAN.



Due to their way of life, ticks can often transmit pathogens from one host to another during a tick bite. There are more types of pathogens involved than in any other parasitic animal group. Also humans are affected by tick-transmitted diseases such as Lyme disease or tick-borne encephalitis (TBE).

The most important vectors in Central Europe are the species of the genus *Ixodes* with the most common native species, the common wood tick (*Ixodes ricinus*).¹⁻³

In recent years, the risk of tick bites has continued to increase, particularly in the summer months, and so-called endemic zones have been established in some areas of Europe.

The effects of an infection with the tick-borne encephalitis virus, range from mild flu-like symptoms to severe neurological deficits and inflammation of the brain and spinal cord.^{1,2}

Serological tests provide information on the antibody status (IgG and IgM) against the TBE virus and can help to confirm an acute infection.³

Serology

Significance

IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

RELIABLE DETECTION OF TBE VIRUS INFECTIONS.

Product Name	TBE/FSME IgG plus ELISA*	TBE/FSME Virus IgM ELISA*
Cat-.Nr.	30114077	30114078
Method	ELISA	ELISA
Kit Size	12x8	12x8
Incubation time	1x1h, 1x30min, 1x15min	1x1h, 1x30min, 1x15min
Standard range	0-300 NTU/mL	Cut-off
Specimen	10 µL serum, plasma (citrate, heparin)	10 µL serum, plasma (citrate, heparin)
Regulatory status	CE IVD	CE IVD

Your benefits:

- Reliable supportive diagnosis of an acute infection
- Easy to automate due to compatible reagents
- Interpretation of results regarding the immune status after vaccination
- Fast sample preparation and processing time
- Excellent performance data

*distributed by Tecan, IBL International GmbH.

Literature:

- 1.Beaute' J., Spiteri G., Warns-Petit E., Zeller H. Tick-borne encephalitis in Europe, 2012 to 2016. Surveillance and outbreak report. Eurosurveillance. Article submitted on 18 Apr 2018 / accepted on 16 Aug 2018 / published on 08 Nov 2018. 2018
- 2.Augustin Estrada-Pena, Frans Jongjean: Ticks feeding on humans: a review of records on human-biting Ixodoidea with special reference to pathogen transmission. In: Experimental and Applied Acarology. September 1999, Band 23, Nr. 9, S. 685-715, PMID 10581710.
- 3.Holzmann, H. (2003). Diagnosis of tick-borne encephalitis. Vaccine, 21, S36-S40.

Manufactured and Distributed by

IBL International GmbH Phone: +49 (0)40-53 28 91-0
 Flughafenstrasse 52a Fax: +49 (0)40-53 28 91-11
 22335 Hamburg Email: IBL@Tecan.com
 Germany www.tecan.com/ibl

Tecan Group Ltd. makes every effort to include accurate and up-to-date information within this publication, however, it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. All mentioned trademarks are protected by law. In general, the trademarks and designs referenced herein are trademarks, or registered trademarks, of Tecan Group Ltd., Männedorf, Switzerland. A complete list may be found at <http://www.tecan.com/trademarks>. Product names and company names that are not contained in the list but are noted herein may be the trademarks of their respective owners. For technical details and detailed procedures of the specifications provided in this document please contact your Tecan representative.

Tecan is in major countries a registered trademark of Tecan Group Ltd., Männedorf, Switzerland.

© 2024 Tecan Trading AG, Switzerland, all rights reserved.

www.tecan.com

