



CONCENTRATED LC-MS KITS.

Consistent performance while reducing costs: With our concentrated LC-MS/MS kits, both laboratory users and purchasers benefit from an innovative solution that combines a sustainable approach with cost-effectiveness, without compromising performance.

Our concentrated kits are specifically designed for professional use in LC-MS/MS analysis. Choose between Starter Kits for initial setup or Solvent Kits for reordering. All components are also available individually, so you can tailor your orders precisely to your needs. The concentrated format allows for straightforward and rapid sample preparation. You save time and can focus on what matters most: precise analyses and reliable results. Our kits guarantee stable and secure storage of reagents. The concentrated solutions minimize the risk of contamination and ensure reproducible results.

By using concentrated kits, you reduce shipping costs thanks to lower weight and volume. Fewer hazardous materials simplify logistics, and compact packaging saves storage space and costs.

Kit Composition

Starter Kit

This Kit includes everything necessary for the first installation of the LC-MS/MS kit on an MS system. Includes all solvents for HPLC and sample extraction, Calibrators, Controls, as well as one HPLC column.

Solvent Kit

The Solvent Kit includes a set of solvents to perform the sample extraction and the HPLC. Calibrators, Controls, Consumables, or an HPLC-Column are not included.

Furthermore, all components can be ordered separately.

Immunosuppressants		Solvent Kit	Starter Kit
Cyclosporin A, Everolimus, Sirolimus, Tacrolimus		30261557	30261558
Deproteinization Solution Concentrated	30257634	1 x 23 ml	1 x 23 ml
Mobile Phase IA Concentrated	30257630	1 x 25 ml	1 x 25 ml
Mobile Phase IB Concentrated	30257631	1 x 25 ml	1 x 25 ml
Mobile Phase II Concentrated	30257633	1 x 12.5 ml	1 x 12.5 ml
Zinc Sulphate	30257635	3 x 25 ml	3 x 25 ml
Tec-Trace Calibrators	30246484	-/-	6 x 2 x 250 µl
Tec-Trace Controls	30246489	-/-	3 x 3 x 250 µl
LC Column Immunosuppressants	30255153	-/-	1 Piece

Disclaimer: For Research use only. Not for use in diagnostic procedures. Distributed by Tecan, IBL International GmbH

Vitamin A and E		Solvent Kit	Starter Kit
		30261564	30261565
Mobile Phase I Concentrated	30256521	1 x 25 mL	1 x 25 mL
Mobile Phase II Concentrated	30256522	1 x 25 mL	1 x 25 mL
Deproteinization solution with internal standard Concentrated	30256520	1 x 20 mL	1 x 20 mL
Tec-Trace Calibrator	30232632	-/-	6 x 2 x 500 µL
Tec-Trace Controls	30232633	-/-	3 x 3 x 500 µL
LC Column Vitamin A & E	30243719	-/-	1 pcs.

Vitamin B1 and B6		Solvent Kit	Starter Kit
		30261534	30261535
Deproteinization Solution Concentrated	30255318	1 x 23.0 mL	1 x 23.0 mL
Mobile Phase I Concentrated	30255314	1 x 25.0 mL	1 x 25.0 mL
Mobile Phase II Concentrated	30255317	1 x 12.5 mL	1 x 12.5 mL
Internal Standard D3	30255319	1 x 6.00 mL	1 x 6.00 mL
Tec-Trace Calibrators	30232634	-/-	6 x 2 x 250 µL
Tec-Trace Controls	30232635	-/-	3 x 3 x 250 µL
LC column Vitamin B1 and B6	30243723	-/-	1 pcs

Methylmalonic Acid		Solvent Kit	Starter Kit
		30261568	30261569
Mobile Phase I Concentrated	30256530	1 x 25 mL	1 x 25 mL
Mobile Phase II Concentrated	30256531	1 x 25 mL	1 x 25 mL
Deproteinization solution	30256529	1 x 16 mL	1 x 16 mL
Solvent	30256533	1 x 20 mL	1 x 20 mL
Internal Standard	30261626	1 x 6 mL	1 x 6 mL
Tec-Trace Calibrator	30232630	-/-	6 x 2 x 500 µL
Tec-Trace Controls	30232631	-/-	3 x 3 x 500 µL
LC Column Methylmalonic Acid	30243761	-/-	1 pcs

Cortisol and Cortisone		Solvent Kit	Starter Kit
		30261570	30261571
Mobile Phase I Concentrated	30256554	1 x 25 mL	1 x 25 mL
Mobile Phase II Concentrated	30256555	1 x 25 mL	1 x 25 mL
Conditioning Solution Concentrated	30256557	1 x 25 mL	1 x 25 mL
Elution Solution Concentrated	30256558	1 x 25 mL	1 x 25 mL
Reconstitution Solution Concentrated	30256569	1 x 10 mL	1 x 10 mL
Washing Solution Concentrated	30256556	2 x 27 mL	2 x 27 mL
Internal standard	30256570	1 x 6 mL	1 x 6 mL
Tec-Trace Calibrator	30232626	-/-	6 x 2 x 1000 µL
Tec-Trace Controls	30232627	-/-	3 x 3 x 1000 µL
LC Column Cortisol and Cortisone	30243736	-/-	1 pcs.

Tecan is in major countries a registered trademark of Tecan Group Ltd., Männedorf, Switzerland.

© 2025 Tecan Trading AG, Switzerland, all rights reserved.

www.tecan.com

