

## ELISA

External Control Material (Urine)



Bio-Rad Lyphochek® Quantitative Urine Control

 USA, Bio-Rad Laboratories Inc.

Update 14.02.2025

IBL Product	Kit Catalog-No.	Bio-Rad Catalog-No.	Bio-Rad Lot	Expiry date	Dilution	Unit	Level 1			Level 2		
							Mean	Lower limit	Upper limit	Mean	Lower limit	Upper limit
<b>5-HIAA ELISA</b>	RE59131	376 / 377 / 375X	63520	2025-05-31	-	mg/L	5.6	3.4	7.9	30.5	24.8	36.2
			63530	2026-01-31	-	mg/L	6.8	5.0	8.6	41.1	22.7	59.4
<b>Adrenalin ELISA</b>	RE59395	376 / 377 / 375X	63520	2025-05-31	-	ng/mL	10.7	7.3	14.1	79.3	35.1	123.4
	RE59242		63530	2026-01-31	-	ng/mL	9.6	7.9	11.3	61.4	34.6	88.1
<b>Dopamine ELISA</b>	RE59395	376 / 377 / 375X	63520	2025-05-31	-	ng/mL	39.9	13.2	66.6	297	220	374
	RE59161		63530	2026-01-31	-	ng/mL	42.2	26.4	58.0	273	211	335
<b>Metanephrine ELISA</b>	RE59181	376 / 377 / 375X	63520	2025-05-31	-	µg/L	90	42	137	326	260	391
			63530	2026-01-31	-	µg/L	129	87	172	423	318	527
			63540	2026-10-31	-	µg/L	113	60	165	386	245	528
<b>Noradrenalin ELISA</b>	RE59395	376 / 377 / 375X	63520	2025-05-31	-	ng/mL	20.2	12.6	27.7	82	67	96
	RE59242		63530	2026-01-31	-	ng/mL	25.0	16.5	33.5	87	60	115
	RE59261											
<b>Normetanephrine ELISA</b>	RE59171	376 / 377 / 375X	63520	2025-05-31	-	µg/L	319	243	396	1175	879	1470
			63530	2026-01-31	-	µg/L	318	197	439	1156	925	1387
			63540	2026-10-31	-	µg/L	275	104	445	1049	622	1477