

Quality Assessment Scheme for Salivary Steroids

DHEA Sample Number 23-21


Your result: The results from all participants can be found in the outstanding table.

Method		LUM	ELISA		MS	
Supplier	all	IBL	IBL	Demeditec	DRG	In-House
Lab. Code	pg/mL	pg/mL	pg/mL	pg/mL	pg/mL	pg/mL
3013	79,1	79,1				
3040	83,6		83,6			
3047	47,8					47,8
3048	115,3		115,3			
3060	112,0		112,0			
3061	102,4		102,4			
6005			[1795]			
13004	121,3		121,3			
15033	89,4					89,4
22236	63,0					
30001	98,0	98,0				
40001	129,0		129,0			
40002	94,3	94,3				
40003-1	101,9					101,9
40013	131,4				131,4	
40015	115,3			115,3		
40016	135,3			135,3		
40024	114,0		114,0			
40027	142,6		142,6			
40029	133,7			133,7		
40036	148,6		148,6			
40041	136,3			136,3		
40049	109,0		109,0			
40055	109,2		109,2			
60006	104,9					104,9

Statistical Results

Method		LUM	ELISA		MS	
Supplier	all	IBL	IBL	Demeditec	DRG	In-House
Mean	109,1	90,5	117,0	130,1	131,4	86,0
SD	24,8	10,0	18,2	10,0		26,3
Mean - 2SD	59,5	70,5	80,6	110,2		33,4
Mean + 2SD	158,7	110,4	153,4	150,1		138,6
Amount	24	3	11	4	1	4

*The graphic is for information purposes only and not validated. Please read your results from the table.



Quality Assessment Scheme for Salivary Steroids

DHEA Sample Number 23-22


Your result: The results from all participants can be found in the outstanding table.

Method		LUM	ELISA		MS	
Supplier	all	IBL	IBL	Demeditec	DRG	In-House
Lab. Code	pg/mL	pg/mL	pg/mL	pg/mL	pg/mL	pg/mL
3013	499,9	499,9				
3040	449,6		449,6			
3047	388,8					388,8
3048	605,2		605,2			
3060	504,7		504,7			
3061	602,4		602,4			
6005			[2940]			
13004	567,8		567,8			
15033	612,7					612,7
22236	570,0					
30001	651,3	651,3				
40001	668,0		668,0			
40002	558,1	558,1				
40003-1	607,3					607,3
40013	750,7				750,7	
40015	1150,0			1150,0		
40016	801,7			801,7		
40024	584,0		584,0			
40027	660,2		660,2			
40029	817,7			817,7		
40036	554,6		554,6			
40041	862,1			862,1		
40049	477,0		477,0			
40055	619,5		619,5			
60006	844,4					844,4

Statistical Results

Method		LUM	ELISA		MS	
Supplier	all	IBL	IBL	Demeditec	DRG	In-House
Mean	642,0	569,8	572,1	907,9	750,7	613,3
SD	164,0	76,4	70,9	163,4		186,1
Mean - 2SD	314,0	417,0	430,3	581,0		241,2
Mean + 2SD	969,9	722,5	713,9	1234,7		985,4
Amount	24	3	11	4	1	4

*The graphic is for information purposes only and not validated. Please read your results from the table.



Quality Assessment Scheme for Salivary Steroids

DHEA Sample Number 23-23

Your result: The results from all participants can be found in the outstanding table.

Method	LUM		ELISA		MS	
Supplier	all	IBL	IBL	Demeditec	DRG	In-House
Lab. Code	pg/mL	pg/mL	pg/mL	pg/mL	pg/mL	pg/mL
3013	320,5	320,5				
3040	389,0		389,0			
3047	279,3					279,3
3048	403,5		403,5			
3060	378,1		378,1			
3061	392,9		392,9			
6005			[1588]			
13004	448,1		448,1			
15033	343,4					343,4
22236	394,0					
30001	420,8	420,8				
40001	480,0		480,0			
40002	358,6	358,6				
40003-1	432,2					432,2
40013	507,2				507,2	
40015	729,2			729,2		
40016	588,2			588,2		
40024	407,0		407,0			
40027	468,7		468,7			
40029	642,8			642,8		
40036	449,1		449,1			
40041	681,3			681,3		
40049	385,0		385,0			
40055	430,2		430,2			
60006	536,0					536,0

Statistical Results

Method	LUM		ELISA		MS	
Supplier	all	IBL	IBL	Demeditec	DRG	In-House
Mean	452,7	366,6	421,1	660,4	507,2	397,7
SD	112,3	50,6	35,8	59,7		111,5
Mean - 2SD	228,1	265,4	349,5	541,0		174,7
Mean + 2SD	677,4	467,8	492,7	779,8		620,7
Amount	24	3	11	4	1	4

*The graphic is for information purposes only and not validated. Please read your results from the table.

