

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Aspergillus fumigatus IgG**

Prod.No.: **30113468**

Lot: **ASPG-061**

Expiry: **2026-11-30**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0680-061	2025-05	2026-11-30
<b>CONJ</b>	1	KG0680-061	2025-05	2026-11-30
<b>CONTROL -</b>	1	NG0680-061	2025-05	2026-11-30
<b>CUT OFF</b>	1	CG0680-061	2025-05	2026-11-30
<b>CONTROL +</b>	1	PG0680-061	2025-05	2026-11-30
<b>WASH BUF 20x</b>	1	W-1112	2025-03	2027-03-31
<b>SOLN STOP</b>	1	S-1043	2025-03	2027-03-31
<b>SUB TMB</b>	1	T-077	2022-09	2027-03-31
<b>DIL</b>	1	D-491	2025-02	2027-02-28

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.014	< 0.200 < Cut-off
<b>CUT OFF</b>	0.648	0.150-1.300
<b>CONTROL +</b>	1.409	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	20.86	11.47 - 30.24

The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-06-17

Approved by:   
 Quality Control

Approved on: 2025-06-18



Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany



IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany