

**GOLD  
STANDARD  
DIAGNOSTICS**



# Certificate of Analysis (CoA)

Product: **Aspergillus fumigatus IgG**

Prod.No.: **30113468**

Lot: **ASPG-064**

Expiry: **2027-07-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0680-064	2026-01	2027-07-31
<b>CONJ</b>	1	KG0680-064	2026-01	2027-07-31
<b>CONTROL -</b>	1	NG0680-064	2026-01	2027-07-31
<b>CUT OFF</b>	1	CG0680-064	2026-01	2027-07-31
<b>CONTROL +</b>	1	PG0680-064	2026-01	2027-07-31
<b>WASH BUF 20x</b>	1	W-1134	2026-02	2028-02-29
<b>SOLN STOP</b>	1	S-1066	2025-12	2027-12-31
<b>SUB TMB</b>	1	T-079	2025-12	2028-12-10
<b>DIL</b>	1	D-499	2026-01	2028-01-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.008	< 0.200 < Cut-off
<b>CUT OFF</b>	0.438	0.150-1.300
<b>CONTROL +</b>	0.900	> Cut-off


Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	23.51	12.93 - 34.08


The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-02-17

Approved by:   
 Quality Control

Approved on: 2026-02-17

 Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany

 IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany