

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Aspergillus fumigatus IgG**

Prod.No.: **30113468**

Lot: **ASPG-067**

Expiry: **2027-12-31**

IFU: **Rev. 01**

| Component           | Number | Lot        | Prod. Date | Expiry Date |
|---------------------|--------|------------|------------|-------------|
| <b>MTP</b>          | 1      | MG0680-067 | 2026-06    | 2027-12-31  |
| <b>CONJ</b>         | 1      | KG0680-067 | 2026-06    | 2027-12-31  |
| <b>CONTROL -</b>    | 1      | NG0680-067 | 2026-06    | 2027-12-31  |
| <b>CUT OFF</b>      | 1      | CG0680-067 | 2026-06    | 2027-12-31  |
| <b>CONTROL +</b>    | 1      | PG0680-067 | 2026-06    | 2027-12-31  |
| <b>WASH BUF 20x</b> | 1      | W-1143     | 2026-05    | 2028-05-31  |
| <b>SOLN STOP</b>    | 1      | S-1075     | 2026-04    | 2028-04-30  |
| <b>SUB TMB</b>      | 1      | T-080      | 2026-02    | 2029-02-23  |
| <b>DIL</b>          | 1      | D-502      | 2026-04    | 2028-04-30  |


| Component        | OD Value | OD Validation Criteria |
|------------------|----------|------------------------|
| <b>CONTROL -</b> | 0.057    | < 0.200 < Cut-off      |
| <b>CUT OFF</b>   | 0.471    | 0.150-1.300            |
| <b>CONTROL +</b> | 0.930    | > Cut-off              |

| Component        | Target Value [NTU] | Target Range [NTU] |
|------------------|--------------------|--------------------|
| <b>CONTROL +</b> | 21.80              | 11.99 - 31.60      |

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-06-15

Approved by:   
Quality Control

Approved on: 2026-06-15



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany