


# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Bordetella pertussis IgM**

Prod.No.: **30113474**

Lot: **BOPM-218**

Expiry: **2027-09-30**

IFU: **Rev. 01**

| Component           | Number | Lot        | Prod. Date | Expiry Date |
|---------------------|--------|------------|------------|-------------|
| <b>MTP</b>          | 1      | MM0030-218 | 2026-03    | 2027-09-30  |
| <b>CONJ</b>         | 1      | KM0030-218 | 2026-03    | 2027-09-30  |
| <b>CONTROL -</b>    | 1      | NM0030-218 | 2026-03    | 2027-09-30  |
| <b>CUT OFF</b>      | 1      | CM0030-218 | 2026-03    | 2027-09-30  |
| <b>CONTROL +</b>    | 1      | PM0030-218 | 2026-03    | 2027-09-30  |
| <b>WASH BUF 20x</b> | 1      | W-1137     | 2026-02    | 2028-02-29  |
| <b>SOLN STOP</b>    | 1      | S-1069     | 2026-02    | 2028-02-29  |
| <b>SUB TMB</b>      | 1      | T-079      | 2025-12    | 2028-12-10  |
| <b>DIL M</b>        | 1      | DM-352     | 2026-03    | 2028-02-29  |

| Component        | OD Value | OD Validation Criteria |
|------------------|----------|------------------------|
| <b>CONTROL -</b> | 0.014    | < 0.200 < Cut-off      |
| <b>CUT OFF</b>   | 0.694    | 0.150-1.300            |
| <b>CONTROL +</b> | 1.085    | > Cut-off              |

| Component        | Target Value [NTU] | Target Range [NTU] |
|------------------|--------------------|--------------------|
| <b>CONTROL +</b> | 17.81              | >10 - 25.82        |

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-03-31

Approved by:   
Quality Control

Approved on: 2026-03-31



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany