


Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Bordetella pertussis toxin (PT) IgA**

Prod.No.: **30113475**

Lot: **BPTA-095**

Expiry: **2027-12-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MA0610-095	2026-06	2027-12-31
CONJ	1	KA0610-095	2026-06	2027-12-31
CAL A	1	SAA0610-095	2026-06	2027-12-31
CAL B	1	SBA0610-095	2026-06	2027-12-31
CAL C	1	SCA0610-095	2026-06	2027-12-31
CAL D	1	SDA0610-095	2026-06	2027-12-31
CONTROL L	1	LCA0610-095	2026-06	2027-12-31
CONTROL H	1	HCA0610-095	2026-06	2027-12-31
WASH BUF 20x	1	W-1144	2026-05	2028-05-31
SOLN STOP	1	S-1076	2026-04	2028-04-30
SUB TMB	1	T-080	2026-02	2029-02-23
DIL	1	D-502	2026-04	2028-04-30

Component	Concentration	OD Value	OD Validation Criteria
CAL A	0 IU/mL	0.014	< 0.200
CAL B	10 IU/mL	0.676	> Standard A
CAL C	25 IU/mL	1.615	> Standard B
CAL D	50 IU/mL	2.458	> 1.000 & > Standard C

Component	Target Value (IU/mL)	Target Range (IU/mL)
CONTROL L	4.14	0 - 10
CONTROL H	30.49	16.77 - 44.20

Final Report issued: 2026-07-06

Approved by: 
Quality Control

Approved on: 2026-07-06



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany