

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Candida albicans IgA**

Prod.No.: **30113617**

Lot: **CANA-202N**

Expiry: **2026-06-30**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MA0060-202	2024-12	2026-06-30
<b>CONJ</b>	1	KA0060-202	2024-12	2026-06-30
<b>CONTROL -</b>	1	NA0060-202	2024-12	2026-06-30
<b>CUT OFF</b>	1	CA0060-202	2024-12	2026-06-30
<b>CONTROL +</b>	1	PA0060-202	2024-12	2026-06-30
<b>WASH BUF 20x</b>	1	W-1101	2024-10	2026-10-31
<b>SOLN STOP</b>	1	S-1032	2024-09	2026-09-30
<b>SUB TMB</b>	1	T-076	2022-05	2026-11-30
<b>DIL</b>	1	D-488	2024-08	2026-08-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.011	< 0.200 < Cut-off
<b>CUT OFF</b>	0.582	0.150-1.300
<b>CONTROL +</b>	1.028	> Cut-off

Final Report issued: 2024-12-17

Approved by:   
Quality Control

Approved on: 2024-12-17



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany