

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Candida albicans IgA**

Prod.No.: **30113617**

Lot: **CANA-212**

Expiry: **2027-10-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MA0060-212	2026-04	2027-10-31
CONJ	1	KA0060-212	2026-04	2027-10-31
CONTROL -	1	NA0060-212	2026-04	2027-10-31
CUT OFF	1	CA0060-212	2026-04	2027-10-31
CONTROL +	1	PA0060-212	2026-04	2027-10-31
WASH BUF 20x	1	W-1139	2026-03	2028-03-31
SOLN STOP	1	S-1071	2026-03	2028-03-31
SUB TMB	1	T-079	2025-12	2028-12-10
DIL	1	D-501	2026-03	2028-03-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.012	< 0.200 < Cut-off
CUT OFF	0.675	0.150-1.300
CONTROL +	1.302	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	18.49	10.17 - 26.81

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-05-04

Approved by: Hoezalp
Quality Control

Approved on: 2026-05-04



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany