

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Echinococcus IgG**

Prod.No.: **30113623**

Lot: **ECHG-096**

Expiry: **2027-05-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0130-096	2025-11	2027-05-31
CONJ	1	KG0130-096	2025-11	2027-05-31
CONTROL -	1	NG0130-096	2025-11	2027-05-31
CUT OFF	1	CG0130-096	2025-11	2027-05-31
CONTROL +	1	PG0130-096	2025-11	2027-05-31
WASH BUF 20x	1	W-1128	2025-11	2027-11-30
SOLN STOP	1	S-1061	2025-11	2027-11-30
SUB TMB	1	T-078	2025-04	2028-04-28
DIL	1	D-496	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
CONTROL -	0.013	< 0.200 < Cut-off
CUT OFF	0.433	0.150-1.300
CONTROL +	0.855	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	22.24	12.23 - 32.24

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-12-15

Approved by: 
Quality Control

Approved on: 2025-12-15



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany