

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Influenza Virus A IgG**

Prod.No.: **30114050**

Lot: **INFAG-153N**

Expiry: **2026-06-30**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0290-153	2024-12	2026-06-30
CONJ	1	KG0290-153	2024-12	2026-06-30
CONTROL -	1	NG0290-153	2024-12	2026-06-30
CUT OFF	1	CG0290-153	2024-12	2026-06-30
CONTROL +	1	PG0290-153	2024-12	2026-06-30
WASH BUF 20x	1	W-1101	2024-10	2026-10-31
SOLN STOP	1	S-1032	2024-09	2026-09-30
SUB TMB	1	T-077	2022-09	2027-03-31
DIL	1	D-488	2024-08	2026-08-31

Component	OD Value	OD Validation Criteria
CONTROL -	0.063	< 0.200 < Cut-off
CUT OFF	0.620	0.150-1.300
CONTROL +	1.251	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	21.24	11.68 - 30.79

The target value and target range of **CONTROL +** may be used optionally.
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-01-14

Approved by: 
 Quality Control

Approved on: 2025-01-14



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany