

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Influenza Virus B IgG**

Prod.No.: **30114053**

Lot: **INFBG-208N**

Expiry: **2026-05-31**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0300-208	2024-11	2026-05-31
<b>CONJ</b>	1	KG0300-208	2024-11	2026-05-31
<b>CONTROL -</b>	1	NG0300-208	2024-11	2026-05-31
<b>CUT OFF</b>	1	CG0300-208	2024-11	2026-05-31
<b>CONTROL +</b>	1	PG0300-208	2024-11	2026-05-31
<b>WASH BUF 20x</b>	1	W-1099	2024-09	2026-09-30
<b>SOLN STOP</b>	1	S-1030	2024-09	2026-09-30
<b>SUB TMB</b>	1	T-076	2022-05	2026-11-30
<b>DIL</b>	1	D-487	2024-07	2026-07-31

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.021	< 0.200 < Cut-off
<b>CUT OFF</b>	0.850	0.150-1.300
<b>CONTROL +</b>	1.408	> Cut-off

Final Report issued: 2024-11-28

Approved by:   
Quality Control

Approved on: 2024-11-28



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany