

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Influenza Virus B IgM**

Prod.No.: **30114054**

Lot: **INFBM-212**

Expiry: **2026-10-31**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MM0300-212	2025-04	2026-10-31
<b>CONJ</b>	1	KM0300-212	2025-04	2026-10-31
<b>CONTROL -</b>	1	NM0300-212	2025-04	2026-10-31
<b>CUT OFF</b>	1	CM0300-212	2025-04	2026-10-31
<b>CONTROL +</b>	1	PM0300-212	2025-04	2026-10-31
<b>WASH BUF 20x</b>	1	W-1109	2025-02	2027-02-28
<b>SOLN STOP</b>	1	S-1040	2025-02	2027-02-28
<b>SUB TMB</b>	1	T-077	2022-09	2027-03-31
<b>DIL M</b>	1	DM-344	2025-04	2027-04-30

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.006	< 0.200 < Cut-off
<b>CUT OFF</b>	0.576	0.150-1.300
<b>CONTROL +</b>	1.071	> Cut-off

Final Report issued: 2025-04-23

Approved by: *S. Herrenschmidt*  
Quality Control

Approved on: 2025-04-23



Gold Standard Diagnostics Frankfurt GmbH  
Waldstrasse 23 A6  
63128 Dietzenbach, Germany



IBL International GmbH  
Flughafenstrasse 52a, 22335 Hamburg  
Germany