

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Mumps Virus IgM**

Prod.No.: **30114066**

Lot: **MUMM-142N**

Expiry: **2026-08-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MM0340-142	2025-02	2026-08-31
CONJ	1	KM0340-142	2025-02	2026-08-31
CONTROL -	1	NM0340-142	2025-02	2026-08-31
CUT OFF	1	CM0340-142	2025-02	2026-08-31
CONTROL +	1	PM0340-142	2025-02	2026-08-31
WASH BUF 20x	1	W-1106	2025-01	2027-01-31
SOLN STOP	1	S-1037	2024-11	2026-11-30
SUB TMB	1	T-077	2022-09	2027-03-31
DIL M	1	DM-342	2024-11	2026-11-30

Component	OD Value	OD Validation Criteria
CONTROL -	0.000	< 0.200 < Cut-off
CUT OFF	0.641	0.150-1.300
CONTROL +	1.331	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	21.59	11.88 - 31.31

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2025-03-11

Approved by: ACr
Quality Control

Approved on: 2025-03-11



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany