

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Clostridium tetani toxin 5S IgG plus**

Prod.No.: **30114074**

Lot: **PTETG-147N**

Expiry: **2026-05-31**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG5043-147	2024-11	2026-05-31
CONJ	1	KG5043-147	2024-11	2026-05-31
CAL A	1	SAG5043-147	2024-11	2026-05-31
CAL B	1	SBG5043-147	2024-11	2026-05-31
CAL C	1	SCG5043-147	2024-11	2026-05-31
CAL D	1	SDG5043-147	2024-11	2026-05-31
CAL E	1	SEG5043-147	2024-11	2026-05-31
CONTROL	1	CG5043-147	2024-11	2026-05-31
WASH BUF 20x	1	W-1100	2024-10	2026-10-31
SOLN STOP	1	S-1031	2024-09	2026-09-30
SUB TMB	1	T-076	2022-05	2026-11-30
DIL	1	D-488	2024-08	2026-08-31

Component	Concentration	OD Value	OD Validation Criteria
CAL A	0.0 IU/mL	0.006	< 0.200
CAL B	0.1 IU/mL	0.289	> 0.050 & > Standard A
CAL C	0.5 IU/mL	0.672	> Standard B
CAL D	1.0 IU/mL	1.046	> Standard C
CAL E	5.0 IU/mL	1.655	> 0.800 & > Standard D

Component	Target Value (IU/mL)	Target Range (IU/mL)
CONTROL	2.13	1.17 - 3.08

Final Report issued: 2024-12-05

Approved by: 

 Quality Control



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany

Approved on: 2024-12-10



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany