

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Clostridium tetani toxin 5S IgG plus**

Prod.No.: **30114074**

Lot: **PTETG-148N**

Expiry: **2026-09-30**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG5043-148	2025-03	2026-09-30
CONJ	1	KG5043-148	2025-03	2026-09-30
CAL A	1	SAG5043-148	2025-03	2026-09-30
CAL B	1	SBG5043-148	2025-03	2026-09-30
CAL C	1	SCG5043-148	2025-03	2026-09-30
CAL D	1	SDG5043-148	2025-03	2026-09-30
CAL E	1	SEG5043-148	2025-03	2026-09-30
CONTROL	1	CG5043-148	2025-03	2026-09-30
WASH BUF 20x	1	W-1107	2025-01	2027-01-31
SOLN STOP	1	S-1038	2024-12	2026-12-31
SUB TMB	1	T-077	2022-09	2027-03-31
DIL	1	D-490	2024-12	2026-12-31

Component	Concentration	OD Value	OD Validation Criteria
CAL A	0.0 IU/mL	0.023	< 0.200
CAL B	0.1 IU/mL	0.423	> 0.050 & > Standard A
CAL C	0.5 IU/mL	0.934	> Standard B
CAL D	1.0 IU/mL	1.377	> Standard C
CAL E	5.0 IU/mL	2.091	> 0.800 & > Standard D

Component	Target Value (IU/mL)	Target Range (IU/mL)
CONTROL	2.31	1.27 - 3.35

Final Report issued: 2025-03-31

Approved by: 
Quality Control



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany

Approved on: 2025-03-31



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany