

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Clostridium tetani toxin 5S IgG plus**
 Prod.No.: **30114074**
 Lot: **PTETG-155**
 Expiry: **2027-09-30**
 IFU: **Rev. 02**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG5043-155	2026-03	2027-09-30
CONJ	1	KG5043-155	2026-03	2027-09-30
CAL A	1	SAG5043-155	2026-03	2027-09-30
CAL B	1	SBG5043-155	2026-03	2027-09-30
CAL C	1	SCG5043-155	2026-03	2027-09-30
CAL D	1	SDG5043-155	2026-03	2027-09-30
CAL E	1	SEG5043-155	2026-03	2027-09-30
CONTROL	1	CG5043-155	2026-03	2027-09-30
WASH BUF 20x	1	W-1139	2026-03	2028-03-31
SOLN STOP	1	S-1071	2026-03	2028-03-31
SUB TMB	1	T-079	2025-12	2028-12-10
DIL	1	D-500	2026-02	2028-02-29

Component	Concentration	OD Value	OD Validation Criteria
CAL A	0.0 IU/mL	0.000	< 0.200
CAL B	0.1 IU/mL	0.297	> 0.050 & > Standard A
CAL C	0.5 IU/mL	0.724	> Standard B
CAL D	1.0 IU/mL	1.114	> Standard C
CAL E	5.0 IU/mL	1.932	> 0.800 & > Standard D

Component	Target Value (IU/mL)	Target Range (IU/mL)
CONTROL	2.08	1.14 - 3.02

Final Report issued: 2026-04-21

Approved by: Hoezalp
 Quality Control

Approved on: 2026-04-21



Gold Standard Diagnostics Frankfurt GmbH
 Waldstrasse 23 A6
 63128 Dietzenbach, Germany



IBL International GmbH
 Flughafenstrasse 52a, 22335 Hamburg
 Germany