

Quality Control Certificate



17beta-Estradiol Saliva ELISA

REF 30121045

LOT EES142


2028-07-31
YYYY-MM-DD

Final Release Results:

	Concentration pg/mL	manual			e.g. automated/ EVOlyzer
		Measured value OD	Range* OD	OD/ODmax	Measured value OD
CAL A	0	2.259	> 1.000	100%	2.778
CAL B	2	1.970	-	87%	2.573
CAL C	4	1.794	-	79%	2.297
CAL D	8	1.407	-	62%	1.782
CAL E	24	0.692	-	31%	0.913
CAL F	64	0.324	-	14%	0.347

*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

Validation criteria (manual/automated)				
	Concentration	Acceptable range		
	found pg/mL	lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	3.1	1.9	3.2	4.5
CONTROL 2	15.8	10.8	15.4	20.0



Kit Composition:

REF		LOT	
30126911	MTP	K913	2028-07
30126912	ENZCONJ	K917	2028-07
30126913	CAL A-F	K914	2028-07
30126914	CONTROL 1	K915	2028-07
30126915	CONTROL 2	K916	2028-07
KEZZ991	TMB SUBS	K610-6	2029-08
KEZZ921	TMB STOP	K919	2029-09
KIZZ01	CONC WASHBUF	K918	2030-03

REF	LOT	
-----	-----	--

Valid edition of the instructions for use: 2023-02

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-10-15

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl