

# Quality Control Certificate




## Reverse T3 (rT3) ELISA

<b>REF</b> 30123153	<b>LOT</b> 253870	<b>2026-11-06</b> YYYY-MM-DD
---------------------	-------------------	---------------------------------

**Final Release Results:**

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.733	> 1.000	100%
CAL B	0.02	2.517	-	92%
CAL C	0.1	1.983	-	73%
CAL D	0.4	1.128	-	41%
CAL E	1	0.635	-	23%
CAL F	2	0.380	-	14%



	Result	Acceptable Range
OD Std. F/% BINDING	0.380/14%	< 20%

	Concentration found ng/mL	Target Value ng/mL	Acceptable range lower limit ng/mL	upper limit ng/mL
Control Low	CS3756	0.17	0.10	0.22
Control High	CS3756	0.44	0.33	0.55

**This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.**

Verified by: Dr. Spangenberg

Validation date: 2025-12-10

**Accepted**

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)