

Quality Control Certificate



Adiponectin ELISA

REF 30126762


LOT 060725

2027-04-17
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	2	0.077	-	3%
CAL B	10	0.372	-	13%
CAL C	30	1.078	-	37%
CAL D	70	2.191	-	75%
CAL E	100	2.906	-	100%

	Target value ng/mL	Acceptable range ng/mL
CONTROL 1	6410	4487-8333
CONTROL 2	13990	9793-18187



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-08-21

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl