


Quality Control Certificate



Transglutaminase IgG ELISA

REF 30132293

LOT 25120107

 2027-05-31
YYYY-MM-DD

Final Release Results:

	Conc.	Measured value	Acceptable range		
	U/mL		lower limit	target value	upper limit
		OD	OD	OD	OD
Calibrator A	0	0.056			< 0.100
Calibrator F	100	2.209	> 1.000		
CONTROL -	-	0.171			≤ 0.300
CONTROL +	-	1.534	≥ 0.700		≤ 2.200

	Validation criteria (manual)			
	Measured value	Acceptable range		
		lower limit	target value	upper limit
	U/mL	U/mL	U/mL	U/mL
CONTROL -	0.9			≤ 2.5
CONTROL +	27	17	27	37

Factor for qualitative assay evaluation	0.31
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2026-02-10

This document has been produced electronically and is valid without a signature.