


Quality Control Certificate



Transglutaminase IgA ELISA

REF 30132295

LOT 26010090

 2027-05-31
YYYY-MM-DD

Final Release Results:

	Conc.	Measured value	Acceptable range		
	U/mL	OD	lower limit OD	target value OD	upper limit OD
Calibrator A	0	0.052			< 0.100
Calibrator F	100	2.362	> 1.000		
CONTROL -	-	0.086			≤ 0.300
CONTROL +	-	1.221	≥ 0.500		≤ 1.600

	Validation criteria (manual)			
	Measured value U/mL	lower limit U/mL	target value U/mL	upper limit U/mL
CONTROL -	0.5			≤ 2.5
CONTROL +	25	16	25	34

Factor for qualitative assay evaluation	0.21
---	------

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Accepted

Validation date: 2026-02-10

This document has been produced electronically and is valid without a signature.