

Quality Control Certificate



HMGB1 Express ELISA

REF 30164033


LOT EHMx134

2026-10-31
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL 1	0	0.035	-	3%
CAL 2	2.5	0.061	-	4%
CAL 3	5.0	0.094	-	7%
CAL 4	10	0.173	-	13%
CAL 5	20	0.332	-	24%
CAL 6	40	0.676	-	49%
CAL 7	80	1.382	-	100%

	Concentration found ng/mL	Acceptable range		
		lower limit ng/mL	target value ng/mL	upper limit ng/mL
CONTROL +	21	14	20	27



Kit Composition:

REF	LOT	Icon	REF	LOT	Icon
30164035 MTP	K290	2026-10	30164042 ENZCONJDIL	K295	2026-10
30164041 LYO CONC ENZCONJ	K291	2026-10	30170004 CONC WASHBUF	K299	2027-01
30164037 CAL	K292	2026-10	30169998 COLREA A	K297	2027-01
30164038 CONTROL +	K293	2026-10	30170002 COLREA B	K298	2027-01
30164040 DILBUF	K294	2026-10	30170003 COLOUR STOP	K300	2027-07

Valid edition of the instructions for use: 2023-07

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2025-09-17

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl