

Quality Control Certificate




Testosterone Saliva Luminescence Immunoassay

REF 30191568 & 30191569	LOT LTSV113	2026-07-31 YYYY-MM-DD
--------------------------------	--------------------	---------------------------------

Final Release Results:

	Concentration pg/mL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	3583028	-	100%
CAL B	6.4	2572193	-	72%
CAL C	16	2244133	-	63%
CAL D	40	1477668	-	41%
CAL E	100	865948	-	24%
CAL F	250	448513	-	13%
CAL G	760	173545	-	5%

	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target pg/mL	upper limit pg/mL
CONTROL 1	26.8	16.9	26.0	35.1
CONTROL 2	182	145	182	218



Kit Composition:

REF		LOT	
KLTE651/9	CONTROL 1	L731	2026-07
KLTE661/9	CONTROL 2	L732	2026-07
KLZZ511	MTP	L728	2027-05
KLTE541/9	CONC ENZCONJ	L859	2026-08
KLTE611/9	CAL A	L730	2026-07
KLTE621/9	CAL B-G	L730	2026-07

REF		LOT	
KLTE521/9	ANTISERUM	L785	2026-07
KLZZ681/9	ASSAYBUF	L721	2026-11
KLZZ601/9	LUMINREAG AP	L813-2	2026-12
KLZZ879	CONC WASHBUF	L727	2026-11

Valid edition of the instructions for use: 2023-12

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: M.Dries

Accepted

Validation date: 2025-02-19

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl