

Quality Control Certificate




sInterleukin-2-Receptor ELISA

REF 30201813**LOT** EIS131 **2027-01-31**
YYYY-MM-DD**Final Release Results:**

<u>Final Release Results:</u>					e.g. automated* DSX
	Concentration U/mL	Measured value OD	Range OD	OD/ODmax	Measured value OD
CAL A	0	0.008	-	0%	0.007
CAL B	7.0	0.167	-	7%	0.163
CAL C	20.0	0.431	-	19%	0.442
CAL D	40.0	0.819	-	35%	0.743
CAL E	80.0	1.472	-	64%	1.399
CAL F	160.0	2.316	-	100%	2.286

*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

	Concentration found U/mL	Acceptable range		
		lower limit U/mL	target U/mL	upper limit U/mL
CONTROL 1	17	10	15	20
CONTROL 2	63	42	63	84

**Kit Composition:**

REF		LOT	
30201806	MTP	K718	2027-11
30201807	CAL A-F	K853	2028-01
30201808	CONTROL 1	K854	2028-01
30201809	CONTROL 2	K855	2028-01
30201810	BIOTIN AB	K597	2027-09
30201811	SAMPLEDIL	K594	2028-02
30201812	ENZCONJ	K859	2028-01
K1ZZZ01	CONC WASHBUF	K596	2029-08
KEZZ991	TMB SUBS	K503-3	2029-04
K1ZZZ71	TMB STOP	K698	2029-03

Valid edition of the instructions for use:

2025-02

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-07-17

Accepted

This document has been produced electronically and is valid without a signature.

**IBL International GmbH**
Flughafenstr. 52A, 22335 Hamburg, GermanyTel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl