

# Quality Control Certificate



## Haemophilus influenzae B IgG ELISA

<b>REF</b> 30217458	<b>LOT</b> IHIB142U	<b>2027-04-30</b> YYYY-MM-DD
---------------------	---------------------	---------------------------------

### Final Release Results:

	Concentration µg/mL	manual		
		Measured value OD	Acceptable range*	OD/ODmax
CAL 1	0.09	0.096	< 0.500	3%
CAL 2	0.17	0.145		5%
CAL 3	0.60	0.534		19%
CAL 4	1.80	1.475		53%
CAL 5	5.00	2.762	> 1.000	100%

\*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions.  
In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

	Concentration found µg/mL	Validation criteria (manual/automated)		
		lower limit µg/mL	target µg/mL	upper limit µg/mL
CONTROL LL	0.53	0.35	0.50	0.65
CONTROL HL	1.19	0.86	1.23	1.59

### Kit Composition:

REF		LOT		REF		LOT	
KIHIB11	MTP	L991	2029-04	KIZZZ71	TMB STOP	L969	2028-04
KIHIB41	ENZCONJ CONC	L964	2027-04				
KIHIB51	CAL 1-5	L965	2028-04				
KIHIB61	CONTROL LL	L966	2028-04				
KIHIB71	CONTROL HL	L967	2028-04				
KIHIB81	DILBUF	L974	2028-04				
KIZZZ01	CONC WASHBUF	L968	2028-10				
KEZZ991	TMB SUBS	L919	2028-06				

Valid edition of the instructions for use: 2022-06

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-09-18

# Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891-0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)