

# Certificate of Analysis (CoA)



**GOLD  
STANDARD  
DIAGNOSTICS**



Product: **Epstein-Barr Virus (EA-D) IgG**  
 Prod.No.: **30219300**  
 Lot: **EAG-006-1**  
 Expiry: **2027-01-31**  
 IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
<b>MTP</b>	1	MG0080-006-1	2025-10	2027-01-31
<b>CONJ</b>	1	KG0080-006-1	2025-10	2027-01-31
<b>CONTROL -</b>	1	NG0080-006-1	2025-10	2027-01-31
<b>CUT OFF</b>	1	CG0080-006-1	2025-10	2027-01-31
<b>CONTROL +</b>	1	PG0080-006-1	2025-10	2027-01-31
<b>WASH BUF 20x</b>	1	W-1124	2025-10	2027-10-31
<b>SOLN STOP</b>	1	S-1055	2025-09	2027-09-30
<b>SUB TMB</b>	1	T-078	2025-04	2028-04-28
<b>DIL</b>	1	D-495	2025-09	2027-09-30

Component	OD Value	OD Validation Criteria
<b>CONTROL -</b>	0.057	< 0.200 < Cut-off
<b>CUT OFF</b>	0.569	0.150-1.300
<b>CONTROL +</b>	1.253	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
<b>CONTROL +</b>	27.55	15.15 - 39.95

The target value and target range of **CONTROL +** may be used optionally.  
 The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-05-29

Approved by:   
 Quality Control

Approved on: 2026-05-29



Gold Standard Diagnostics Frankfurt GmbH  
 Waldstrasse 23 A6  
 63128 Dietzenbach, Germany



IBL International GmbH  
 Flughafenstrasse 52a, 22335 Hamburg  
 Germany