

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Yersinia enterocolitica IgA**

Prod.No.: **30219306**

Lot: **YERA-023**

Expiry: **2027-04-30**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MA0990-023	2026-07	2027-07-31
CONJ	1	KA0990-023	2026-07	2027-07-31
CONTROL -	1	NA0990-023	2026-07	2027-07-31
CUT OFF	1	CA0990-023	2026-07	2027-07-31
CONTROL +	1	PA0990-023	2026-07	2027-04-30
WASH BUF 20x	1	W-1146	2026-06	2028-06-30
SOLN STOP	1	S-1078	2026-06	2028-06-30
SUB TMB	1	T-080	2026-02	2029-02-23
DIL	1	D-503	2026-06	2028-06-30


Component	OD Value	OD Validation Criteria
CONTROL -	0.011	< 0.200 < Cut-off
CUT OFF	0.654	0.150-1.300
CONTROL +	1.222	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	19.15	10.53 - 27.76

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-08-10

Approved by: 
Quality Control

Approved on: 2026-08-10



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany