

Certificate of Analysis (CoA)



**GOLD
STANDARD
DIAGNOSTICS**



Product: **Yersinia enterocolitica IgG**

Prod.No.: **30219308**

Lot: **YERG-020**

Expiry: **2027-07-31**

IFU: **Rev. 01**

Component	Number	Lot	Prod. Date	Expiry Date
MTP	1	MG0990-020	2026-04	2027-07-31
CONJ	1	KG0990-020	2026-04	2027-07-31
CONTROL -	1	NG0990-020	2026-04	2027-07-31
CUT OFF	1	CG0990-020	2026-04	2027-07-31
CONTROL +	1	PG0990-020	2026-04	2027-07-31
WASH BUF 20x	1	W-1139	2026-03	2028-03-31
SOLN STOP	1	S-1073	2026-04	2028-04-30
SUB TMB	1	T-080	2026-02	2029-02-23
DIL	1	D-500	2026-02	2028-02-29

Component	OD Value	OD Validation Criteria
CONTROL -	0.012	< 0.200 < Cut-off
CUT OFF	0.459	0.150-1.300
CONTROL +	0.799	> Cut-off

Component	Target Value [NTU]	Target Range [NTU]
CONTROL +	17.52	>10.00 - 25.40

The target value and target range of **CONTROL +** may be used optionally.

The positive control is valid if the absorbance value of **CONTROL +** > **CUT OFF**.

Final Report issued: 2026-05-04

Approved by: 
Quality Control

Approved on: 2026-05-04



Gold Standard Diagnostics Frankfurt GmbH
Waldstrasse 23 A6
63128 Dietzenbach, Germany



IBL International GmbH
Flughafenstrasse 52a, 22335 Hamburg
Germany