

Quality Control Certificate




Cortisol Saliva Luminescence Immunoassay

REF 30221150 & 30221151	LOT LCL206	2026-06-30 YYYY-MM-DD
--------------------------------	-------------------	---------------------------------

Final Release Results:

	Concentration µg/dL	Measured value RLU	Range RLU	RLU/RLUmax
CAL A	0	27854133	-	100%
CAL B	0.015	19849585	-	71%
CAL C	0.03	16640385	-	60%
CAL D	0.10	8494630	-	30%
CAL E	0.40	2747158	-	10%
CAL F	1.00	1069198	-	4%
CAL G	3.20	270898	-	1%

	Concentration found µg/dL	Acceptable range		
		lower limit µg/dL	target µg/dL	upper limit µg/dL
CONTROL 1	0.082	0.052	0.081	0.109
CONTROL 2	0.76	0.51	0.73	0.94



Kit Composition:

REF	LOT		REF	LOT			
KLCL651/9	CONTROL 1	M187D-5D	2026-12	KLZZ641/9	LUMINREAG 1	L753	2026-10
KLCL661/9	CONTROL 2	M188D-5D	2026-12	KLZZ651/9	LUMINREAG 2	L754	2026-10
KLCL511	MTP	L715	2026-12				
KLCL541/9	ENZCONJ	L773	2026-12				
KLCL611/9	CAL A	M185-5D	2026-12				
KLCL621/9	CAL B-G	M186-5D	2026-12				
KLZZ709	CONC WASHBUF	L711	2026-12				

Valid edition of the instructions for use: 2022-10

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2024-12-17

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl