

# Quality Control Certificate




## Angiotensin I (PRA) ELISA

<b>REF</b> DB52011	<b>LOT</b> 261240	<b>2027-03-24</b> YYYY-MM-DD
--------------------	-------------------	---------------------------------

**Final Release Results:**

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	0	2.748	> 1.200	100%
CAL B	0.2	2.522	-	92%
CAL C	0.5	2.257	-	82%
CAL D	1.5	1.694	-	62%
CAL E	4	1.008	-	37%
CAL F	10	0.536	-	20%
CAL G	25	0.266	-	10%
CAL H	60	0.155	-	6%



	Result	Acceptable Range
OD Std. H/% BINDING	0.155/6%	< 10%

	Concentration found ng/mL	Target Value ng/mL	Acceptable range lower limit ng/mL	upper limit ng/mL
KIT CONTROL Low CS3891	0.93	0.93	0.70	1.16
KIT CONTROL High CS3891	9.74	9.88	7.41	12.35

**This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.**

Verified by: Dr. Spangenberg

Validation date: **Accepted**

2026-05-08

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)