

# Quality Control Certificate



## Pregnenolone ELISA

REF DB52031


LOT 261930

2027-11-30  
YYYY-MM-DD

### Final Release Results:

|       | Concentration<br>ng/mL | Measured<br>value<br>OD | Range<br>OD | OD/ODmax |
|-------|------------------------|-------------------------|-------------|----------|
| CAL A | 0                      | 2.315                   | > 1.500     | 100%     |
| CAL B | 0.1                    | 1.671                   | -           | 72%      |
| CAL C | 0.4                    | 1.239                   | -           | 54%      |
| CAL D | 1.6                    | 0.719                   | -           | 31%      |
| CAL E | 6.4                    | 0.392                   | -           | 17%      |
| CAL F | 25.6                   | 0.197                   | -           | 9%       |

|                     | Result   | Acceptable Range |
|---------------------|----------|------------------|
| OD Std. F/% BINDING | 0.197/9% | < 30%            |



|              | Concentration<br>found<br>ng/mL | Target<br>Value<br>ng/mL | Acceptable range<br>lower limit<br>ng/mL | upper limit<br>ng/mL |
|--------------|---------------------------------|--------------------------|------------------------------------------|----------------------|
| Control Low  | CS3937                          | 1.17                     | 0.87                                     | 1.45                 |
| Control High | CS3937                          | 4.14                     | 3.37                                     | 5.61                 |

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-06-18

**Accepted**

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)