

Quality Control Certificate




Estrone ELISA

| | | |
|--------------------|-------------------|---------------------------------|
| REF DB52051 | LOT 253990 | 2027-05-19 YYYY-MM-DD |
|--------------------|-------------------|---------------------------------|

Final Release Results:

| | Concentration pg/mL | Measured value OD | Range OD | OD/ODmax |
|-------|------------------------|-------------------------|-------------|----------|
| CAL A | 0 | 2.368 | > 1.100 | 100% |
| CAL B | 20 | 2.080 | - | 88% |
| CAL C | 60 | 1.718 | - | 73% |
| CAL D | 200 | 1.206 | - | 51% |
| CAL E | 600 | 0.708 | - | 30% |
| CAL F | 2000 | 0.363 | - | 15% |



| | Result | Acceptable Range |
|-----------------------------|--------|------------------|
| Std. F % BINDING (OD/ODmax) | 15% | < 20% |

| | Concentration found pg/mL | Target Value pg/mL | Acceptable range lower limit pg/mL | upper limit pg/mL |
|-------------------------|---------------------------------|--------------------------|--|----------------------|
| KIT CONTROL Low CS3766 | 82.1 | 85.8 | 51.5 | 120.1 |
| KIT CONTROL High CS3766 | 485.7 | 485.9 | 364.4 | 607.4 |

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-12-15

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl