

# Quality Control Certificate




## Androstane-diol- Glucuronide ELISA

|            |                |            |               |                                 |
|------------|----------------|------------|---------------|---------------------------------|
| <b>REF</b> | <b>DB52171</b> | <b>LOT</b> | <b>252760</b> | <b>2026-12-31</b><br>YYYY-MM-DD |
|------------|----------------|------------|---------------|---------------------------------|

### Final Release Results:

|       | Concentration<br>ng/mL | Measured<br>value<br>OD | Range<br>OD | OD/ODmax |
|-------|------------------------|-------------------------|-------------|----------|
| CAL A | 0                      | 2.701                   | > 2.000     | 100%     |
| CAL B | 0.25                   | 1.816                   | -           | 67%      |
| CAL C | 1                      | 1.152                   | -           | 43%      |
| CAL D | 3                      | 0.648                   | -           | 24%      |
| CAL E | 10                     | 0.284                   | -           | 11%      |
| CAL F | 50                     | 0.093                   | -           | 3%       |

|                         | Result   | Acceptable Range |
|-------------------------|----------|------------------|
| OD Standard F/% Binding | 0.093/3% | < 12%            |



|                         | Concentration<br>found<br>ng/mL | Target<br>Value<br>ng/mL | Acceptable range<br>lower limit<br>ng/mL | upper limit<br>ng/mL |
|-------------------------|---------------------------------|--------------------------|--|----------------------|
| KIT CONTROL High CS3611 | 8.02                            | 8.33                     | 6.25                                     | 10.41                |
| KIT CONTROL Low CS3611  | 1.87                            | 1.88                     | 1.41                                     | 2.35                 |

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-09-09

# Accepted

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)