

Quality Control Certificate



Free Testosteron(e) ELISA

REF DB52181

LOT EFT146


2028-06-30
YYYY-MM-DD

Final Release Results:

	Concentration pg/mL	manual			e.g. automated/ DSX
		Measured value OD	Range*	OD/ODmax	Measured value OD
CAL A	0	2.260	-	100%	2.106
CAL B	0.25	2.118	-	94%	1.903
CAL C	2.5	1.615	-	71%	1.449
CAL D	10	1.011	-	45%	0.904
CAL E	25	0.572	-	25%	0.506
CAL F	100	0.146	-	6%	0.135

*During automated procedure the optical densities of calibrators and controls may rise due to temperature and specific device conditions. In case of none conformance with the OD ranges stated on the QC Certificate please contact IBL-International or use the validation criteria only!

Validation criteria (manual/automated)				
	Concentration found pg/mL	Acceptable range		
		lower limit pg/mL	target value pg/mL	upper limit pg/mL
CONTROL 1	2.5	1.5	2.5	3.5
CONTROL 2	17.7	12.5	17.8	23.2



Kit Composition:

REF		LOT	
KEFT511	MTP	K835	2028-06
KEFT541	ENZCONJ	K839	2028-06
KEFT611	CAL A-F	K836	2028-06
KEFT651	CONTROL 1	K837	2028-06
KEFT661	CONTROL 2	K838	2028-06
KIZZ01	CONC WASHBUF	K840	2029-12
KEZZ991	TMB SUBS	K503-4	2029-04
KEZZ921	TMB STOP	K841	2029-06

REF	LOT	
-----	-----	--

Valid edition of the instructions for use: 2025-07

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: T. Borchert

Validation date: 2026-08-18

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl