

Quality Control Certificate



IgA Saliva ELISA

REF DM59171


LOT 6543A

2027-04-30
YYYY-MM-DD

Final Release Results:

	Calibrator concentration µg/mL	Measured value OD	Range OD	OD/ODmax
CAL 0	0	0.099	-	3%
CAL 1	6.9	0.278	> 0.110	9%
CAL 2	62	1.110	> 0.445	36%
CAL 3	132	1.765	> 0.705	57%
CAL 4	400	3.103	> 1.240	100%

	Concentration found µg/mL	Acceptable range		
		lower limit µg/mL	target µg/mL	upper limit µg/mL
CONTROL	99.2	55.4	85.2	115.0




This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-02-13

Accepted

This document has been produced electronically and is valid without a signature.

 IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl