

Quality Control Certificate




IgA Saliva ELISA

REF DM59171	LOT 6578A	2027-06-30 YYYY-MM-DD
--------------------	------------------	---------------------------------

Final Release Results:

	Calibrator concentration µg/mL	Measured value OD	Range OD	OD/ODmax
CAL 0	0	0.095	-	4%
CAL 1	6.9	0.237	> 0.095	10%
CAL 2	62	0.972	> 0.390	41%
CAL 3	132	1.568	> 0.625	67%
CAL 4	400	2.348	> 0.940	100%



	Concentration found µg/mL	Acceptable range		
		lower limit µg/mL	target µg/mL	upper limit µg/mL
CONTROL	74.6	53.1	81.7	110.3

This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-04-08

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl