

Quality Control Certificate



Leptin ELISA

REF MD53001


LOT 030624-1

2026-07-02
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	1	0.113	-	4%
CAL B	10	0.926	-	31%
CAL C	25	1.638	-	55%
CAL D	50	2.254	-	76%
CAL E	100	2.961	-	100%

	Concentration found ng/mL	Acceptable range ng/mL
CONTROL 1	11.94	8.36 - 15.52
CONTROL 2	24.68	17.28 - 32.08



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2024-08-20

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl