

# Quality Control Certificate



## Leptin ELISA

REF MD53001


LOT 031125-2

2027-12-02  
YYYY-MM-DD

### Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	1	0.092	-	3%
CAL B	10	0.839	-	30%
CAL C	25	1.508	-	55%
CAL D	50	2.161	-	78%
CAL E	100	2.764	-	100%

	Concentration found ng/mL	Acceptable range ng/mL
CONTROL 1	11.00	7.70 - 14.30
CONTROL 2	24.53	17.17 - 31.89



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-02-04

# Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)