

# Quality Control Certificate




## Leptin ELISA

<b>REF</b> MD53001	<b>LOT</b> 130426-2	<b>2028-05-12</b> YYYY-MM-DD
--------------------	---------------------	---------------------------------

### Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	1	0.103	-	4%
CAL B	10	0.872	-	30%
CAL C	25	1.604	-	56%
CAL D	50	2.307	-	80%
CAL E	100	2.885	-	100%

	Concentration found ng/mL	Acceptable range ng/mL
CONTROL 1	11.82	8.27 - 15.37
CONTROL 2	26.20	18.34 - 34.06



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-06-04

# Accepted

This document has been produced electronically and is valid without a signature.



**IBL International GmbH**  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)