

Quality Control Certificate



Leptin ELISA

REF MD53001


LOT 200525-2

2027-06-19
YYYY-MM-DD

Final Release Results:

| | Concentration ng/mL | Measured value OD | Range OD | OD/ODmax |
|-------|------------------------|-------------------------|-------------|----------|
| CAL A | 1 | 0.099 | - | 4% |
| CAL B | 10 | 0.797 | - | 29% |
| CAL C | 25 | 1.514 | - | 55% |
| CAL D | 50 | 2.050 | - | 74% |
| CAL E | 100 | 2.752 | - | 100% |

| | Concentration found ng/mL | Acceptable range ng/mL |
|-----------|---------------------------------|---------------------------|
| CONTROL 1 | 11.48 | 8.04 - 14.92 |
| CONTROL 2 | 24.42 | 17.09 - 31.75 |



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2025-10-10

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl