

Quality Control Certificate



IGF-1 ELISA

REF MD58011


LOT 070725

2027-08-06
YYYY-MM-DD

Final Release Results:

	Concentration ng/mL	Measured value OD	Range OD	OD/ODmax
CAL A	2	0.075	-	3%
CAL B	5	0.261	-	11%
CAL C	15	0.925	-	39%
CAL D	30	1.692	-	72%
CAL E	50	2.362	-	100%

In-house control samples		
	Target value ng/mL	Acceptable range ng/mL
CONTROL 1	99	80 - 119
CONTROL 2	453	363-544



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-06-24

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.:
E-MAIL:
WEB:

+ 49 (0) 40 532891 -0 Fax: -11
ibl@tecan.com
www.tecan.com/ibl