

Quality Control Certificate



IGFBP-3 ELISA

REF MD58031


LOT 170326

2028-04-16
YYYY-MM-DD

Final Release Results:

	Concentration	Measured value	Range	OD/ODmax
	ng/mL	OD	OD	
CAL A	0.4	0.048	-	2%
CAL B	2	0.219	-	11%
CAL C	6	0.578	-	30%
CAL D	15	1.225	-	64%
CAL E	30	1.929	-	100%

	Target value	Acceptable range
	ng/mL	ng/mL
CONTROL 1	2247	1798-2697
CONTROL 2	4260	3408-5113



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-06-18

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl