

# Quality Control Certificate



## Androstenedione RIA (CT)

REF MG11091

LOT 262104K


2026-07-31  
YYYY-MM-DD

### Final Release Results:

		Measured value	
Total Activity	TA	53892	cpm
Maximum Binding/TA	B0/TA	29.36	%

	Concentration	Measured value	Range	B/B0
	ng/mL	cpm	cpm	
CAL 0	0.0	15824	-	100%
CAL 1	0.11	12852	-	81%
CAL 2	0.31	10569	-	67%
CAL 3	1.04	6773	-	43%
CAL 4	2.84	3852	-	24%
CAL 5	10.0	1528	-	10%

	Measured value	Acceptable range		
		lower limit	target value	upper limit
	ng/mL	ng/mL	ng/mL	ng/mL
CONTROL 1	0.72	0.44	0.65	0.86
CONTROL 2	2.43	1.73	2.47	3.21



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-05-27

**Accepted**

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)