

Quality Control Certificate



17-OH-Progesterone RIA (CT)

REF MG12181

LOT 261904K


2026-07-17
YYYY-MM-DD

Final Release Results:

		Measured value	
Total Activity	TA	76772	cpm
Maximum Binding/TA	B0/TA	32.99	%

	Concentration ng/mL	Measured value cpm	Range cpm	B/B0
CAL 0	0	25324	-	100%
CAL 1	0.22	19139	-	76%
CAL 2	0.49	13962	-	55%
CAL 3	1.53	5738	-	23%
CAL 4	5.03	1730	-	7%
CAL 5	16.3	584	-	2%

	Measured value ng/mL	Acceptable range		
		lower limit ng/mL	target value ng/mL	upper limit ng/mL
CONTROL 1	1.00	0.60	0.94	1.28
CONTROL 2	2.76	1.75	2.73	3.71



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-05-11

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl