

# Quality Control Certificate



## 17-OH-Progesterone RIA (CT)

REF MG12181

LOT 262505K


2026-08-28  
YYYY-MM-DD

### Final Release Results:

|                    |       | Measured value |     |
|--------------------|-------|----------------|-----|
| Total Activity     | TA    | 83439          | cpm |
| Maximum Binding/TA | B0/TA | 35.56          | %   |

|       | Concentration | Measured value | Range | B/B0 |
|-------|---------------|----------------|-------|------|
|       | ng/mL         | cpm            | cpm   |      |
| CAL 0 | 0             | 29667          | -     | 100% |
| CAL 1 | 0.22          | 22361          | -     | 75%  |
| CAL 2 | 0.49          | 15997          | -     | 54%  |
| CAL 3 | 1.53          | 7177           | -     | 24%  |
| CAL 4 | 5.03          | 2139           | -     | 7%   |
| CAL 5 | 16.3          | 712            | -     | 2%   |

|           | Measured value | Acceptable range |              |             |
|-----------|----------------|------------------|--------------|-------------|
|           | ng/mL          | lower limit      | target value | upper limit |
|           | ng/mL          | ng/mL            | ng/mL        | ng/mL       |
| CONTROL 1 | 1.04           | 0.67             | 1.04         | 1.41        |
| CONTROL 2 | 2.61           | 1.61             | 2.52         | 3.43        |



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-06-22

**Accepted**

This document has been produced electronically and is valid without a signature.



IBL International GmbH  
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11  
E-MAIL: [ibl@tecan.com](mailto:ibl@tecan.com)  
WEB: [www.tecan.com/ibl](http://www.tecan.com/ibl)