

Quality Control Certificate



Aldosterone RIA (CT)


REF MG13051	LOT 261705K	2026-06-26 YYYY-MM-DD
--------------------	--------------------	---------------------------------

Final Release Results:

		Measured value	
Total Activity	TA	42442	cpm
Maximum Binding/TA	B0/TA	66.08	%

	Concentration pg/mL	Measured value cpm	Range cpm	B/B0
CAL 0	0.0	28047	-	100%
CAL 1	26.5	23725	-	85%
CAL 2	52.9	20641	-	74%
CAL 3	100.0	17587	-	63%
CAL 4	305.0	11574	-	41%
CAL 5	659.0	7918	-	28%
CAL 6	1542.0	5088	-	18%

	Measured value pg/mL	Acceptable range		
		lower limit pg/mL	target value pg/mL	upper limit pg/mL
CONTROL 1	189.35	135.00	193.00	251.00



This product has been tested successfully by the Quality Control Department and was released for sale according to the existing specifications.

Verified by: Dr. Spangenberg

Validation date: 2026-04-28

Accepted

This document has been produced electronically and is valid without a signature.



IBL International GmbH
Flughafenstr. 52A, 22335 Hamburg, Germany

Tel.: + 49 (0) 40 532891 -0 Fax: -11
E-MAIL: ibl@tecan.com
WEB: www.tecan.com/ibl